



GET BACK TO YOU.

The Spanner® Prostatic Stent has been designed to open the prostatic urethra and allow patients to maintain urine flow and allow volitional voiding by reducing urethral resistance within the prostate.

The Spanner® allows Normal micturition achieved with bladder and sphincter; Improves urination efficiency by decreasing urine left in the bladder; Reduces rate of bacterial colonization compared to the Foley catheter; Maintains sexual activity and is Well tolerated by patients. It also Improves IPSS scores and is preferred by patients over catheterization.

srsmedical.com
999.999.9999

SRS
MEDICAL